

High Content Image Analysis Using a Laser Scanning Microplate Cytometer

Sarah Payne, Ben Schenker Andrew Goulter and Paul Wylie
TTP LabTech Ltd, Melbourn Science Park, Melbourn, Hertfordshire, UK.

Abstract

Microscope-based high-content instruments offer high optical resolution, however, the limited field of view afforded by their objective lenses can mean lengthy read times for some assays, especially where multiple image capture per well is required. Wherever possible, only a small percentage of the total number of cells present in a test well are analysed to keep plate read times at a minimum, which may not always be ideal. Microplate laser-scanning cytometers, such as TTP LabTech's Acumen[®] X3, are an alternative platform. This technology offers higher throughput than imaging based systems processing 300,000 wells of data in 24 hours in 1536 plates. The Acumen[®] X3 is capable of generating cytometric data and TIFF images (8 or 16-bit) simultaneously. The images represent the whole well and correlate with those captured using a 20x microscope objective.


The images produced can be subjected to the full range of algorithms contained in commercially available image analysis software. In most cases where image analysis and cytometric analysis have been validated, both methods were comparable where single cells were present. For a limited number of applications, however, secondary information provided by image analysis may increase understanding. Such assays include angiogenesis, cell colony formation and tissue scanning.

Microplate cytometry offers whole well, high throughput capability with the increased flexibility of simultaneous cytometric and TIFF image outputs and presents a comprehensive solution for high content analysis studies.

Conclusion


- TIFF image export feature extends the application of the Acumen[®] X3 into more complex high content assays requiring image processing.
- Images generated are equivalent to those from a 20X microscope objective.
- Removes the requirement for time consuming generation of multiple images with rare event detection assays.
- There is no need for image stitching steps prior to image analysis on large objects such as cell colonies, tubules, *C.elegans* and tissue sections.
- The speed to capture whole well images is not reduced with higher density plates.
- Exported TIFF files are compatible with many commercial imaging packages.

1 Comparison of High Content Instrumentation Optics




The large field of view (400 mm² (20 x 20 mm)) is far greater than that offered by microscope-based CCD imagers (~1 mm² for a x10 objective). The small field of view on microscope-based CCD imagers is only sufficient to obtain resolved images of around 100 cells at once. The application of laser scanning over a large area means that analysis is performed on an area, not a well basis. This equates to the simultaneous scanning of 4, 16 and 64 wells in 96, 384 and 1536 well format, respectively. Reconfiguration of assays into higher density plate formats results in a concomitant increase in throughput up to 300,000 samples per day in 1536 well microplates.

2 Large field of View Images from an Acumen[®] X3



- A. Whole well image of a viral infectivity assay
B. HUVEC cells stained with calcein-AM
C. Whole section of rat duodenum labeled with a non-specific lipid stain


3 Demonstration of the Resolution of Images from Acumen[®] X3




Mitotic Index: Cells were stained with anti-phospho-histone H3 (FITC) and the plate scanned on an Acumen[®] X3 to generate 8 bit images. The cells were pseudo coloured green. The circular arrangement of DNA within the nucleus is clearly visible upon compound treatment demonstrating that the images generated are approximately equivalent to a 20X microscope objective.

HeLa cells were stained with Hoechst and the plate scanned on an Acumen[®] X3 microplate cytometer. The image shown is an 8 bit TIFF file of PMT measurements. The resolution is sufficient to distinguish daughter cells and dividing cells within the cell population.

4 Flow Chart for Image Analysis of Protein Translocation Assay




5 Image Analysis of a Kinase Nuclear Translocation Assay



The assay was run in a 96 well plate and an image of the whole well area was obtained. The scan time was 25 minutes per plate. The results showed that 25% of cells were found to have translocated kinase. This compared with the data obtained using standard Acumen[®] X3 software with a scan time of 10 minutes per plate. (Data was analysed using Image-Pro Plus version 6.1)

6 Application for Rare Event Detection: Stem Cell Research



Stem cells were scanned using an Acumen[®] X3 microplate cytometer. A & B, cell proliferation using anti-BrdU antibody (Alexa 488 secondary) and PI nuclear counterstain. C & D, stem cell differentiation using calcein-AM (green) and a selective reporter gene (red). Cytometric data were exported as 8-bit TIFF images and image-processed to identify and count differentiated cells within each cluster (Media Cybernetics, Berkshire, UK). Data were supplied by Epistem Ltd, Manchester, UK.



TTP LabTech Ltd., Melbourn Science Park, Melbourn, Hertfordshire, SG8 6EE, UK. Phone +44 1763 262626
Visit us at www.ttplabtech.com/acumen